



# Malaria P.f. /P.v. Antigen Rapid Test Kit (Colloidal Gold Method)

The Malaria P.f. /P.v. Antigen Rapid Test Kit(Colloidal Gold Method) is a rapid chromatographic immunoassay designed for the qualitative detection of P. falciparum (P.f.) and P.vivax (P.v.) antigens in human whole blood. Early diagnosis and treatment of malaria reduces disease, prevents deaths and contributes to reducing transmission. WHO recommends that all suspected cases of malaria be confirmed using parasite-based a rapid diagnostic test. Diagnostic testing enables health providers to swiftly distinguish between malarial and non-malarial fevers, facilitating appropriate treatment.

## Test Procedure

- **Specimen:** Capillary and venous whole blood.
- **Specimen collection:** Specimens of all types are collected by conventional methods.
- **Test method:**



## Specifications

Components	Test card/Sample diluent/ Dropper
Pack Size	25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C~30°C ;Valid for 24 months
Test samples	Human whole blood
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

## Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Type
Malaria P.f. /P.v. Antigen Rapid Test Kit(Colloidal Gold Method)	BSK12S1S	25 Tests/Kit	1 × 3 mL/Tube	Human whole blood
	BSK12M1S	50 Tests/Kit	2 × 3 mL/Tube	