



# Typhoid IgG/IgM Rapid Test Kit (Colloidal Gold Method)

The Typhoid IgG/IgM Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection and differentiation of IgG and IgM Typhoid antibodies produced against Salmonella Typhi (*S. typhi*) in human whole blood, serum or plasma. This easy-to-use format with rapid results allows point-of-care testing, providing information to assist with treatment and patient management decisions.

## Test Procedure

- **Specimen:** Whole blood, serum or plasma.
- **Specimen collection:** Specimens of all types are collected by conventional methods.
- **Test method:**



## Specifications

| Components          | Test card/Sample diluent/ Sample Droppers |
|---------------------|---|
| Pack Size           | 25Tests/Kit; 50 Tests/Kit                 |
| Storage Conditions  | 2°C~30°C ;Valid for 24 months             |
| Test samples        | Whole blood, Serum or Plasma              |
| Reaction conditions | 15~30°C 15~20 minutes                     |
| Special instrument  | No need                                   |

## Product Informations

| Product Name  | Cat.No   | Package      | Sample Diluent | Sample Droppers | Sample Type                    |
|---|----------|--------------|----------------|-----------------|--------------------------------|
| Typhoid IgG/IgM Rapid Test Kit<br>(Colloidal Gold Method) | BSK11S1S | 25 Tests/Kit | 1 × 3 mL/Tube  | 25              | Whole blood,<br>serum & plasma |
|   | BSK11M1S | 50 Tests/Kit | 2 × 3 mL/Tube  | 50              |                                |